

Freezing and Storage of Sperm



**AUSTRALIAN
CONCEPT
FERTILITY
CENTRE**

WHAT DOES THIS INVOLVE?

The freezing of sperm has been performed since the early 1950's and has been undertaken at Concept Fertility Centre since the early 1980's. A semen sample, collected at Concept Fertility Centre is mixed with a cryoprotectant (similar to antifreeze in a car) and drawn into thin plastic straws, each of about ½ cc in volume. Normally between 5 and 20 straws may be prepared from each single ejaculate. The straws are coloured and a coloured plug is used to seal the end. These colour combinations are unique for each man. The man's name, Concept file number and sample number are recorded on each straw. These straws are then placed into a programmable freezing machine and liquid is pumped into it at a controlled rate. This allows the sperm straws to be frozen in a scientifically controlled manner.

At the completion of the freezing programme, the straws of semen are stored in liquid nitrogen at - 196 C in a state of suspended animation. Sperm can stay in this situation for many years with no adverse effect.



WHY FREEZE SPERM?

Prior to having a vasectomy
Prior to chemotherapy or radiation therapy.

When the male partner is going to be away during fertility treatment.

HOW IS SPERM FROZEN?

A semen sample is produced and then mixed with a special freezing media

This mixture is then drawn up into thin, coloured plastic straws (a bit like drinking straws) which are sealed with a coloured plug

The straws are labelled with name and a unique file number and placed into a controlled rate freezer machine

Once frozen the straws are stored at -196°C in liquid nitrogen

WHO WOULD STORE SPERM?

The freezing and storage of sperm is undertaken for many reasons:


Pre-vasectomy/Pre-chemotherapy

Men may store sperm before undergoing vasectomy or chemotherapy as an insurance policy. Vasectomies should be considered as irreversible because repair and reconnection of the vas sometime after vasectomies usually is accompanied by a significant reduction in normal fertilizing capacity of sperm.

Some chemotherapy can result in damage to a man's future ability to produce sperm. Freezing sperm before treatment can be a way to circumvent this potential problem.

Management of Assisted Reproductive Treatment

Some men who are not going to be readily available when their partners are undergoing treatment for assisted conception or men who find semen collection on demand a difficult exercise, could consider the storage of sperm a convenient option.



MANAGEMENT OF SPERM STORAGE

Men who choose to store their sperm need to be aware of the costs. Concept Fertility Centre charges for the freezing and storage of 2 (two) samples. (Please contact Concept Fertility Centre to discuss the costs at the time). Further fees may be incurred after 10 years storage

All men who wish to store sperm need to undergo blood tests for HIV, Hepatitis B and Hepatitis C.

Request forms to store semen must be signed by men before storage, and this consent will last for 5 (five) years up to a maximum of 15 (fifteen) years. Special applications will be required to extend storage beyond 15 (fifteen) years.

SEMEN COLLECTION

Semen samples are collected in the men's room or at home. For convenience it is important to book the men's room at a suitable time (Monday-Friday, 7.30am -2 pm). We advise abstaining from sexual relations for 3 days prior to collection. Because the men's room is used by others it is important not to be late for the appointment.

Concept will inform you if there is any concern from the laboratory about the quality of the sperm or if the sperm is unsuitable for freezing.

It is important to inform Concept of any change of address. If men are not contactable for renewal of the consent and application for continued storage after (15) fifteen years, the samples will be discarded. Under the Human Reproductive Technology Act 1991, no posthumous use of semen is permitted. In the event of death, Concept will need to be informed so that decisions can be made regarding the disposal of the stored semen.

USE OF SPERM

The sperm will be thawed and used only with specific written consent of the man who owns the stored semen.

FURTHER INFORMATION

Should you require any further information please call the Coordinator Doctor at Concept.

Stages involved in sperm freezing

Undergo blood tests for HIV, Hepatitis B & C
Sign a semen storage consent form

Collect a semen sample at Concept or at home

Semen is frozen

How long can semen remain in storage?

The Human Reproductive Technology Act 1991 allows semen to remain in storage for a maximum of 15 years. Consent to continue storage must be given every 5 years. In special circumstances storage extension could be granted.

Change of address?

Please let Concept Fertility Centre know your new contact

Questions and Notes

CONTACT US:

ACFC-KARACHI

F-6/1 Block-8 Near Teen Talwar Behind Gulf Mall Clifton Karachi
Tel: 021-35810049-50 Cell: 0321-2110099

ACFC-LAHORE

15-C Ghalib Road Gulberg II Near Hafeez Centre Lahore
Tel: 042-35880831-34 Cell: 0321-2233074

info@australianconceptfertility.com.pk | info@conceptfertility.com.pk
www.australianconceptfertility.com.pk | www.conceptfertility.com.pk